

QUICK GUIDE

Allow the test components to reach room temperature before use.

1. STEP 1: Sample collection


NOTE! Collection and application of the sample must be swift, in order to avoid sample coagulation in the micropipette.

- Twist off the lancet cap until you feel it separates from the device. This activates the lancet. Do not pull the cap – just twist it
- Press the open end of the lancet on the patient's fingertip and press the release button.
- Wipe away the first drop of blood (using a sterile gauze pad or cotton).
- Collect the blood sample:
 - Hold the blood collection micropipette horizontally, with the air vent free.
 - Touch the blood sample with the tip of the tube. Do not touch or scrape the skin.
 - Fill the micropipette to the fill line, avoiding air bubbles.



2. STEP 2: Test procedure

- Slowly dispense the blood sample (20 μ l) onto the test cassette.
- Wait 30–60 seconds until the blood has been absorbed.
- Add two drops of dilution buffer into the same window. If the dropper is not used, add 100 μ l of dilution buffer.



3. STEP 3: Reading of the results

- Incubate for 10 minutes
- Interpret the results (C = control line, T = test line)
- Do not read the result after 10 min



Negative



Positive

